



Technical Update • February 2015

Cleveland Clinic Laboratories is dedicated to keeping you updated and informed about recent testing changes. That's why we are happy to provide this technical update on a monthly basis.

Recently changed tests are **bolded**, and could include revisions to methodology, reference range, days performed or CPT code. For your convenience, tests are listed alphabetically and the order and billing codes are provided.

If you wish to compare the new information with previous test information, refer to the Test Directory, which can be accessed at clevelandcliniclabs.com.

Deleted tests and new tests are listed separately. Please update your database as necessary. For additional detail, contact Client Services at 216.444.5755 or 800.628.6816 or via email at clientservices@ccf.org.

Test Update Page #	Summary of Changes by Test Name	Order Code	Billing Code	Name Change	New Test	Specimen Special Information	Component Requirement	Methodology Change	Reference Range	Days Performed/Reported	CPT
3	17-Ketosteroids, 24 Hr Urine										
12	Adenovirus Antibody, IgG										
12	Adenovirus Antibody, IgM										
3	AFB Culture Only										
3	Alpha-1-Antitrypsin										
3	Amikacin, Post Dose										
3	Amikacin, Pre Dose										
4	Amikacin, Random										
4	Amitriptyline/Nortriptyline										
4	Amphetamine Confirmation, Urine										
4	Amphetamine Confirmation Quant, Serum/Plasma										
11	Anti-C1q Antibody, IgG										
4	Anti-Streptolysin O										
4	Barbiturates Confirmation, Urine										
4	Benzodiazepines Confirmation, Urine										
4	Blau Syndrome NOD2/CARD15 Complete Gene Analysis										
4	Bromide										
4	Buprenorphine and Metabolite, Confirm/Quant, Serum/Plasma										
4	Catecholamines, Fractionated, Plasma										
12	Catecholamines Fractionated Tilt										
5	Cathartic Laxative, Stool										
5	Ceruloplasmin										
5	Clobazam										
5	Clonazepam & Metabolite, Urine										
5	Date Rape Panel										
5	Disopyramide (Norpace)										
5	Drug Analysis, Comprehensive										
5	Epidermal Antibodies with Reflex to Titer										

Test Update Page #	Summary of Changes by Test Name	Order Code	Billing Code	Name Change	New Test	Specimen Test Discontinued	Special Information	Component Requirement	Component Change	Methodology	Days Performed/Reported Reference Range	CPT
5	Ethosuximide											
5	Ethyl Glucuronide, Urine reflex to Confirm/Quant											
5	FISH Insight Analysis											
5	FTA Antibodies CSF											
5	Fungal Culture and Smear (Non Dermal)											
6	Fungal Culture and Smear Hair, Skin, Nail											
6	Gamma-Hydroxybutyric Acid, Serum											
6	Hepatitis C Virus FibroSURE											
12	Hepatitis C Virus RNA Quantitative bDNA											
12	HIV 1 Quantitative bDNA											
6	HLA-A, B and C											
6	HLA-A29											
6	HLA-B27 PCR											
6	HLA-DR/DQ											
6	HPV Genotypes 16, 18/45											
6	Imipramine/Desipramine											
7	KIT (D816V) Mutation by PCR											
7	LSD, Urine											
7	Meperidine and Normeperidine											
7	Methsuximide/Normethsuximide											
7	Mexiletine											
7	Nocardia Culture Only											
8-9	Organic Acids Urine, Quant											
10	Oxycodone Confirmation, Urine											
12	Paraldehyde & Acetaldehyde											
10	Perphenazine											
10	Phencyclidine Confirmation, Urine											
10	Propafenone											
11	Reticulocyte, Hemoglobin											
10	Sm/RNP Antibody											
10	Sotalol											
10	Streptozyme											
10	Synthetic Cannabinoid Metabolite, Screen with Confirmation											
10	Tapentadol and Metabolite Confirm/Quantitation, Urine											
10	Testosterone, Free, Adult Males by ED/LC-MS/MS											
10	Testosterone, Free/Total, Males by ED/LC-MSMS											
10	Tobramycin, Post Dose											
11	Tobramycin, Pre Dose											
11	Tobramycin, Random											
11	VWF Exon 28 Genotyping*											

*=Test Modification from January Technical Update

Test Changes

Test Name	Order Code	Billing Code	Change	Effective Date
17-Ketosteroids, 24 Hr Urine	U17K	83582	CPT: 83586x1, 82570x1	4/3/2015
AFB Culture Only	AFCO	76091	<p>Specimen Requirements: 10.0 mL Blood - Isolator 10 microbial tube (1.5 mL minimum); Ambient. Recommended volume for adults is 10.0 mL. For pediatric patients, draw 1.5 mL blood into a 1.5 mL Isolator tube.</p> <p>OR - 2.0 mL Bone marrow - Isolator tube (aerobic only - 1.0 mL minimum); Draw 1-2 mL bone marrow using a heparinized syringe.</p> <p>Transfer to a pediatric Isolator tube.</p> <p>Alternate - 2.0 mL Sodium or Lithium heparin (Green - 1.0 mL minimum); Ambient. Acceptable container for bone marrow specimens. Transport promptly to Microbiology.</p> <p>Special Information: An AFB stain will not be performed. Patient preparation: For blood cultures, select vein to use. Wipe off venipuncture site using a 70% alcohol pad. Apply Chloraprep to the skin over the selected venipuncture site and apply using up and down and back and forth strokes for a full 30 seconds. Allow the site to dry completely for 30-60 seconds. Swab septum of Isolator tube or Myco/F bottle with a Chloraprep or 70% alcohol. Draw 10 mL into adult Isolator tube, 1.5 mL into Pediatric Isolator tube or 5 mL if direct draw into the Myco/F bottle. After inoculation, clean septum with alcohol swab. Transport to Microbiology within 4-6 hours is recommended.</p> <p>Clinical Information: An AFB Culture only test should be performed to identify an infection due to mycobacteria in blood or bone marrow specimens. Broth medium will be utilized for culturing blood or bone marrow sites. Identification of positive cultures will be performed utilizing current methodologies. Susceptibility testing will be performed on significant isolates. Additional charges may apply. A single negative culture does not rule the presence of a mycobacterial infection.</p> <p>Stability: Frozen: Unacceptable Refrigerated: Unacceptable Ambient: < 8 Hours</p>	2/16/2015
Alpha-1-Antitrypsin	AAT	30030	<p>Specimen Requirements: 1.0 mL Plasma - lithium heparin (Light Green); Refrigerated.</p> <p>Alternate - 1.0 mL Serum - SST (Gold); Refrigerated.</p> <p>Methodology: Immunoturbidometric Assay</p> <p>Reported: 8 Hours</p> <p>Stability: Frozen: 3 Months Refrigerated: 3 Months Ambient: 7 Days</p>	4/6/2015
Amikacin, Post Dose	AMIKPO	52010	<p>Specimen Requirements: 1.0 mL Plasma - sodium or lithium heparin (Green - 0.5 mL minimum); Refrigerated. Collect 30 minutes after completion of infusion. Centrifuge. Transfer plasma/serum to a clean, tightly sealed tube. Refrigerate. Freeze if storage/transport time will be longer than 24 hours.</p> <p>Alternate - 1.0 mL Serum - SST (Gold - 0.5 mL minimum); Refrigerated.</p> <p>Methodology: Fluorescence Polarization Immunoassay (FPIA)</p> <p>Reported: 8 Hours</p> <p>Stability: Frozen: 1 Year Refrigerated: 48 Hours Ambient: 8 Hours</p>	4/2/2015
Amikacin, Pre Dose	AMIKPR	52008	<p>Specimen Requirements: 1.0 mL Plasma - sodium or lithium heparin (Green - 0.5 mL minimum); Refrigerated. Collect 5-90 minutes before next infusion. Centrifuge. Transfer plasma/serum to a clean, tightly sealed tube. Refrigerate. Freeze if storage/transport time will be longer than 24 hours.</p> <p>Alternate - 1.0 mL Serum - SST (Gold - 0.5 mL minimum); Refrigerated.</p> <p>Methodology: Fluorescence Polarization Immunoassay (FPIA)</p> <p>Reported: 8 Hours</p> <p>Stability: Frozen: 1 Year Refrigerated: 48 Hours Ambient: 8 Hours</p>	4/2/2015

Test Changes (cont.)

Test Name	Order Code	Billing Code	Change	Effective Date
Amikacin, Random	AMIKRA	76075	Specimen Requirements: 1.0 mL Plasma - sodium or lithium heparin (Green - 0.5 mL minimum); Refrigerated. Centrifuge. Transfer plasma/serum to a clean, tightly sealed tube. Refrigerate. Freeze if storage/transport time exceeds 24 hours. Alternate - 1.0 mL Serum - SST (Gold - 0.5 mL minimum); Refrigerated. Methodology: Fluorescence Polarization Immunoassay (FPIA) Reported: 8 Hours Stability: Frozen: 1 Year Refrigerated: 48 Hours Ambient: 8 Hours	4/2/2015
Amitriptyline/Nortriptyline	AMINOR	82138	CPT: 80335x1	1/1/2015
Amphetamine Confirmation, Urine	UAMPC	83603	CPT: 80324x1	1/5/2015
Amphetamine Confirmation Quant, Serum/Plasma	AMPCQ	88679	CPT: 80326x1	1/1/2015
Anti-Streptolysin O	ASO	86060	Methodology: Immunoturbidometric Assay Reported: 8 Hours Stability: Frozen: 6 Months Refrigerated: 2 Days Ambient: 2 Days	4/6/2015
Barbiturates Confirmation, Urine	UBARBC	90310	CPT: 80345x1	1/1/2015
Benzodiazepines Confirmation, Urine	UBENZC	83370	CPT: 80346x1	1/1/2015
Blau Syndrome NOD2/CARD15 Complete Gene Analysis	BLAU	88320	CPT: 81479x1	1/1/2015
Bromide	BROM	82290	Special Information: Draw specimen prior to next dose - at steady state concentration. Please provide the following information if available 1.) Dose - List drug amount and include the units of measure 2.) Route - List the route of administration (IV, oral, etc.) 3.) Dose Frequency - Indicate how often the dose is administered (per day, per week, as needed, etc.) 4.) Type of Draw - Indicate the type of blood draw (Peak, Trough, Random, etc.) Reported: 2-6 Days Component: Test build may need to be modified.	2/17/2015
Buprenorphine and Metabolite, Confirm/Quant, Serum/Plasma	SBUP	89772	CPT: 80348X1	1/1/2015
Catecholamines, Fractionated, Plasma	PLCAT	41000	Specimen Requirements: 4.0 mL Plasma - sodium heparin (Green - 2.5 mL minimum); Frozen. Freeze ASAP. Draw specimen in a pre-chilled green-top Vacutainer. Plasma should be separated within 30 minutes of collection and then frozen immediately at -20 degrees C. Alternate - 4.0 mL Plasma - lithium heparin (Green - 2.5 mL minimum); Frozen. Freeze ASAP. Draw specimen in a pre-chilled green-top Vacutainer. Plasma should be separated within 30 minutes of collection and then frozen immediately at -20 degrees C. Methodology: High Performance Liquid Chromatography with Electrochemical Detection Clinical Information: The three catecholamines (norepinephrine, epinephrine, and dopamine) are the principal secretory products of neural tissue. Clinically, the measurement of circulation catecholamines is valuable in diagnosis of catecholamine secreting tumors associated chiefly with hypertension (pheochromocytomas, neuroblastomas, and gangliomas) and with the evaluation of orthostatic hypotension. Days Performed: Monday-Friday Reported: 4-7 Days Component: Add Total Catecholamines (CATCOL): Total (Norepinephrine+Epinephrine=Dopamine) Supine: 123-671 pg/ml Upright: 242-1125 pg/ml Stability: Frozen: 30 Days Refrigerated: 6 Hours Ambient: 6 Hours	2/10/2015

Test Changes (cont.)

Test Name	Order Code	Billing Code	Change	Effective Date
Cathartic Laxative, Stool	STCATH	83334	CPT: 83735x1, 84100x1, 80375x1	1/1/2015
Ceruloplasmin	CERULO	40009	Specimen Requirements: 1.0 mL Plasma - lithium heparin (Light Green); Frozen. Alternate - 1.0 mL Serum - SST (Gold); Frozen. Methodology: Immunoturbidometric Assay Reported: 8 Hours Stability: Frozen: 4 Weeks Refrigerated: 3 Days Ambient: 8 Hours	4/6/2015
Clobazam	CLOBAZ	90152	CPT: 80339x1	1/1/2015
Clonazepam & Metabolite, Urine	UCLONO	83859	CPT: 80346x1	1/1/2015
Date Rape Panel	UDRPAN	82125	CPT: 80304x1, 80301x1	1/1/2015
Disopyramide (Norpace)	DISOP	34032	Special Information: Draw specimen prior to next dose - at steady state concentration. Please provide the following information if available: 1.) Dose - List drug amount and include the units of measure 2.) Route - List the route of administration (IV, oral, etc.) 3.) Dose Frequency - Indicate how often the dose is administered (per day, per week, as needed, etc.) 4.) Type of Draw - Indicate the type of blood draw (Peak, Trough, Random, etc.) Component: Test build may need to be modified.	2/17/2015
Drug Analysis, Comprehensive	DRANCO	82053	Specimen Requirements: 10.0 mL Whole blood - potassium oxalate/sodium fluoride (Gray); Refrigerated. THIS ASSAY REQUIRES MULTIPLE SPECIMEN TYPES. AND 30.0 mL Urine, random - clean container; Refrigerated. THIS ASSAY REQUIRES MULTIPLE SPECIMEN TYPES. Methodology: Gas Chromatography Mass Spectrometry (GCMS) Reported: 6-11 Days Stability: Frozen: 6 Months Refrigerated: 2 Weeks Ambient: 3 days CPT: 83992x1, 82570x1, 80377x1, 80373x1, 80372x1, 80371x1, 80370x1, 80368x1, 80367x1, 80366x1, 80365x1, 80364x1, 80361x1, 80360x1, 80359x1, 80385x1, 80357x1, 80355x1, 80354x1, 80353x1, 80348x1, 80346x1, 80344x1, 80341x1, 80338x1, 80337x1, 80334x1, 80331x1, 80301x1, 80326x1, 80302x1	2/24/2015
Epidermal Antibodies with Reflex to Titer	EPIABS	89577	Days Performed: Wednesday, Friday Reported: 4-7 Days	4/2/2015
Ethosuximide	ETHOS	82692	Special Information: Draw specimen prior to next dose - at steady state concentration. Please provide the following information if available: 1.) Dose - List drug amount and include the units of measure 2.) Route - List the route of administration (IV, oral, etc.) 3.) Dose Frequency - Indicate how often the dose is administered (per day, per week, as needed, etc.) 4.) Type of Draw - Indicate the type of blood draw (Peak, Trough, Random, etc.) Component: Test build may need to be modified.	2/17/2015
Ethyl Glucuronide, Urine reflex to Confirm/Quant	UEGLUC	89509	CPT: 80321x1	1/1/2015
FISH Insight Analysis	ISIGHT	82933	Specimen Requirements: 20.0 mL Fluid, amniotic - sterile container (18.0 mL minimum); Refrigerate. Do not centrifuge for any reason. Alternate - 1.0-2.0 mL Whole blood - sodium heparin (Green); Ambient. Pubescent or newborn blood.	4/2/2015
FTA Antibodies CSF	FTACSF	79231	Clinical Information: The significance of a reactive FTA-ABS CSF test is unknown. The CSF from persons treated in the secondary or latent stages of syphilis and without signs of neurosyphilis may be reactive. A nonreactive result in the FTA-ABS CSF test suggests the absence of neurosyphilis. Reported: 2-3 Days	2/17/2015
Fungal Culture and Smear (Non Dermal)	FCULSM	77917	CPT: 87210x1 AND 87102X1	4/1/2015

Test Changes (cont.)

Test Name	Order Code	Billing Code	Change	Effective Date
Fungal Culture and Smear Hair, Skin, Nail	FHSNSM	89653	Special Information: Hair, skin, nail and scalp are the only acceptable specimen types. For other sources, please use order code FCUL for culture or FCULSM when requesting both fungal culture and smear. Test includes culture for yeasts and molds. Additional billing is applied for identification and susceptibility testing. CPT codes vary based on methodology. CPT: 87220x1, 87101X1	4/1/2015
Gamma-Hydroxybutyric Acid, Serum	GHBSER	82415	CPT: 80304x1	1/1/2015
Hepatitis C Virus FibroSURE	HCVSUR	84198	Specimen Requirements: 3.0 mL Serum - SST (Gold); Frozen. Centrifuge, aliquot into 2 tubes, one containing 2.5 mL serum and the other 0.5 mL serum, and freeze both tubes ASAP. Alternate - 3.0 mL Serum - no additive (Red); Frozen. Centrifuge, aliquot into 2 tubes, one with 2.5 mL serum and the other with 0.5 mL serum, and freeze both tubes ASAP. Methodology: Nephelometry (NEPH), Colorimetric, Kinetic, Colorimetry, Immunologic Days Performed: Varies, Monday-Friday Reported: 17-21 Days	4/2/2015
HLA-A, B and C	HLABC	82817	Specimen Requirements: 7.0 mL Whole blood - ACD A or B (Yellow); Ambient. Reported: 7-10 Days	4/2/2015
HLA-A29	HLAA29	88634	Specimen Requirements: 5.0 mL Whole blood - EDTA (Lavender); Ambient. Alternate - 7.0 mL Whole blood ACD A or B (Yellow); Ambient. Clinical Information: Birdshot retinochoroidopathy (BSCR) is a rare subtype of idiopathic posterior uveitis with distinct clinical characteristics that can lead to severe visual impairment. BSCR has the strongest documented HLA association for a human disease with > 95% of patients carrying HLA-A29. The frequency of HLA-A29 varies by ethnic group and could be up to 10% in some US populations. Determinations of HLA-A29 is of diagnostic significance in BSCR. Reported: 7-10 Days Primary Name: HLA-A29	4/2/2015
HLA-B27 PCR	B27PCR	83080	Specimen Requirements: 5.0 mL Whole blood - EDTA (Lavender - 5.0 mL minimum); Ambient. Alternate - 7.0 mL Whole blood ACD A or B (Yellow - 5.0 mL minimum); Ambient. Collect one ACD A or B yellow top tube. Methodology: Sequence Specific Oligonucleotide Probe (SSOP) Clinical Information: HLA-B27 is strongly associated with ankylosing spondylitis (AS). HLA-B27 is also associated with other seronegative arthropathies such as Reiter syndrome and psoriatic arthritis as well as extra-articular diseases such as anterior uveitis and inflammatory bowel disease. Greater than 90% of patients with AS are HLA-B27 positive. The frequency of HLA-B27 varies by ethnic group but generally <10% in most US populations. HLA-B27 associated susceptibility to AS varies by population and HLA-B27 alleles detected. Some alleles such as B*27:05 are associated with high AS susceptibility while others such as B*27:06 and B*27:09 are associated with low susceptibility. HLA-B27 allele typing is recommended in HLA-B27 positive cases. Days Performed: Monday-Friday Reported: 7-10 Days Stability: Frozen: Unacceptable Refrigerated: 1 Week Ambient: 1 Week	4/2/2015
HLA-DR/DQ	HLADR	83639	Specimen Requirements: 5.0 mL Whole blood - EDTA (Lavender); Ambient. Alternate - 7.0 mL Whole blood ACD A or B (Yellow); Ambient. Days Performed: Monday-Friday Reported: 7-10 Days Primary Name: HLA-DR/DQ	4/2/2015
HPV Genotypes 16, 18/45	HPVGNO	89544	CPT: 87625x1	1/1/2015
Imipramine/Desipramine	IMIDES	34044	CPT: 80335x1	1/1/2015

Test Changes (cont.)

Test Name	Order Code	Billing Code	Change	Effective Date
KIT (D816V) Mutation by PCR	KITMST	84159	<p>Specimen Requirements: 5.0 mL Whole blood - EDTA (Lavender - 1.0 mL minimum); Refrigerated.</p> <p>Alternate - 3.0 mL Bone marrow - EDTA (Lavender - 1.0 mL minimum); Refrigerated.</p> <ul style="list-style-type: none"> - 100 mg Tissue, fresh - sterile container; Frozen, ASAP. Collect 100 mg or 0.5-2.0 cm of fresh tissue and freeze immediately. - Paraffin block, formalin fixed; Ambient. Submit formalin fixed (10 percent neutral buffered formalin), paraffin embed FFPE tumor tissue. Protect from excessive heat. Transport tissue block or four 10-micron shavings. <p>Methodology: Polymerase Chain Reaction (PCR)</p> <p>Clinical Information: Aids in the diagnosis of mastocytosis. Provide prognostic and predictive information for tyrosine kinase inhibitor (TKI) therapy planning.</p> <p>Days Performed: Monday, Wednesday, Friday</p> <p>Reported: 3-8 Days</p> <p>Component: Test build may need to be modified.</p> <p>Stability: Frozen: Whole blood/Bone marrow: Unacceptable Fresh Tissue: 1 Year FFPE tumor tissue: Unacceptable Refrigerated: Whole blood/Bone marrow: 5 Days Fresh tissue: 2 Hours FFPE tumor tissue: Indefinitely Ambient: Whole blood/Bone marrow: 24 Hours Fresh tissue: Unacceptable FFPE tumor tissue: Indefinitely</p> <p>Primary Name: KIT (D816V) Mutation by PCR</p>	4/2/2015
LSD, Urine	ULSD	88129	CPT: 80302x1	1/1/2015
Meperidine & Normeperidine	MEPNO	80337	CPT: 80362x1	1/1/2015
Methsuximide/ Normethsuximide	METHSU	83627	<p>Special Information: Draw specimen prior to next dose - at steady state concentration. Please provide the following information if available:</p> <ol style="list-style-type: none"> 1.) Dose - List drug amount and include the units of measure 2.) Route - List the route of administration (IV, oral, etc.) 3.) Dose Frequency - Indicate how often the dose is administered (per day, per week, as needed, etc.) 4.) Type of Draw - Indicate the type of blood draw (Peak, Trough, Random, etc.). <p>Component: Test build may need to be modified.</p>	2/17/2015
Mexiletine	MEX	75504	<p>Special Information: Draw specimen prior to next dose - at steady state concentration. Please provide the following information if available:</p> <ol style="list-style-type: none"> 1.) Dose - List drug amount and include the units of measure 2.) Route - List the route of administration (IV, oral, etc.) 3.) Dose Frequency - Indicate how often the dose is administered (per day, per week, as needed, etc.) 4.) Type of Draw - Indicate the type of blood draw (Peak, Trough, Random, etc.). <p>Component: Test build may need to be modified.</p>	2/17/2015
Nocardia Culture Only	NOCARC	77756	<p>Special Information: Indicate specimen source on requisition. Identification of any Nocardia species isolated will be performed, but susceptibility testing will only be performed upon request. Susceptibility testing will be performed by outside reference lab. Additional billing is applied for identification and susceptibility testing. CPT codes vary based on methodology.</p>	4/1/2015

Test Changes (cont.)

Test Name	Order Code	Billing Code	Change	Effective Date
Organic Acids Urine, Quant	UORA	89797	Reference Ranges: 3-Methylglutaconate: 0-29 Days: 0.0-5.2 umol/mmolCr 30-364 Days: 0.0-7.9 umol/mmolCr 1-99 Years: 0.0-2.0 umol/mmolCr ButyrylGlycine, Ur: 0-29 Days: 0.0-2.7 umol/mmolCr 30-364 Days: 0.0-2.7 umol/mmolCr 1-99 Years: 0.0-0.7 umol/mmolCr 3MethylGlutarate, Ur: 0-29 Days: 0.0-0.9 umol/mmolCr 30-364 Days: 0.1-1.4 umol/mmolCr 1-99 Years: 0.0-0.6 umol/mmolCr Glutarate, Urine: 0-29 Days: 0.0-6.2 umol/mmolCr 30-364 Days: 0.4-6.6 umol/mmolCr 1-99 Years: 0.0-1.4 umol/mmolCr Pyruvate, Urine: 0-29 Days: 0.4-5.7 umol/mmolCr 30-364 Days: 0.2-13.2 umol/mmolCr 1-99 Years: 0.1-2.6 umol/mmolCr SuberylGlycine, Ur: 0-29 Days: 0.0 umol/mmolCr 30-364 Days: 0.0 umol/mmolCr 1-99 Years: 0.0 umol/mmolCr N-AcetylTyrosine, Ur: 0-29 Days: 0.0-3.7 umol/mmolCr 30-364 Days: 0.0-4.9 umol/mmolCr 1-99 Years: 0.0-1.2 umol/mmolCr 4OHPhenylLactate, Ur: 0-29 Days: 0.4-112.6 umol/mmolCr 30-364 Days: 1.1-88.2 umol/mmolCr 1-99 Years: 1.3-23.0 umol/mmolCr Sebacic Acid, Urine: 0-29 Days: 0.0-4.6 umol/mmolCr 30-364 Days: 0.0-12.4 umol/mmolCr 1-99 Years: 0.0-0.3 umol/mmolCr MethylCitrate, Urine: 0-29 Days: 0.0-16.8 umol/mmolCr 30-364 Days: 1.7-25.4 umol/mmolCr 1-99 Years: 1.0-13.9 umol/mmolCr IsoCitric Acid, Urine: 0-29 Days: 2.4-297.4 umol/mmolCr 30-364 Days: 43.9-537.3 umol/mmolCr 1-99 Years: 9.1-271.9 umol/mmolCr Aconitate, Urine: 0-29 Days: 5.9-161.3 umol/mmolCr 30-364 Days: 26.3-330.8 umol/mmolCr 1-99 Years: 8.5-109.6 umol/mmolCr 2-OxoAdipic Acid, Ur: 0-29 Days: 0.0-0.0 umol/mmolCr 30-364 Days: 0.0-2.8 umol/mmolCr 1-99 Years: 0.0-3.3 umol/mmolCr SuccinylAcetone, Ur: 0-29 Days: <0.4 umol/mmolCr 30-364 Days: <0.4 umol/mmolCr 1-99 Years: <0.4 umol/mmolCr Suberic Acid, Urine: 0-29 Days: 0.0-23.1 umol/mmolCr 30-364 Days: 2.4-30.7 umol/mmolCr 1-99 Years: 0.0-7.4 umol/mmolCr N-AcetylAsparticAcid: 0-29 Days: 0.2-40.0 umol/mmolCr 30-364 Days: 0.0-69.0 umol/mmolCr 1-99 Years: 0.1-8.9 umol/mmolCr 4OHPhenylAcetate, Ur: 0-29 Days: 3.1-146.6 umol/mmolCr 30-364 Days: 14.3-569.5 umol/mmolCr 1-99 Years: 5.7-147.5 umol/mmolCr HexanoylGlycine, Ur: 0-29 Days: 0.0-0.5 umol/mmolCr 30-364 Days: 0.0-0.3 umol/mmolCr 1-99 Years: 0.0-0.1 umol/mmolCr a-KetoGlutarate, Ur: 0-29 Days: 0.0-403.2 umol/mmolCr 30-364 Days: 0.2-355.8 umol/mmolCr 1-99 Years: 0.2-42.7 umol/mmolCr 2HydroxyGlutaricAcid: 0-29 Days: 3.0-49.1 umol/mmolCr 30-364 Days: 4.4-66.9 umol/mmolCr 1-99 Years: 0.6-17.7 umol/mmolCr 3HydroxyGlutaricAcid: 0-29 Days: 0.0-0.7 umol/mmolCr 30-364 Days: 0.0-2.5 umol/mmolCr 1-99 Years: 0.0-0.7 umol/mmolCr	4/2/2015

Test Changes (cont.)

Test Name	Order Code	Billing Code	Change	Effective Date
Organic Acids Urine, Quant continued			3MECrotonylGlycine, Ur: 0-29 Days: <0.3 umol/mmolCr	4/2/2015
			30-364 Days: <0.3 umol/mmolCr	
			1-99 Years: <0.3 umol/mmolCr	
5-Oxo-Proline, Urine:			0-29 Days: 0.0-7.3 umol/mmolCr	
			30-364 Days: 0.0-7.8 umol/mmolCr	
			1-99 Years: 0.4-3.1 umol/mmolCr	
Adipic Acid, Urine:			0-29 Days: 0.5-52.5 umol/mmolCr	
			30-364 Days: 5.0-53.8 umol/mmolCr	
			1-99 Years: 0.3-9.2 umol/mmolCr	
Malate, Urine:			0-29 Days: 0.4-12.8 umol/mmolCr	
			30-364 Days: 0.9-11.1 umol/mmolCr	
			1-99 Years: 0.0-1.1 umol/mmolCr	
2MEButyrylGlycine, Ur:			0-29 Days: 0.0-1.0 umol/mmolCr	
			30-364 Days: 0.0-0.9 umol/mmolCr	
			1-99 Years: 0.0-0.4 umol/mmolCr	
IsoButyrylGlycine, Ur:			0-29 Days: 0.0-1.1 umol/mmolCr	
			30-364 Days: 0.0-1.1 umol/mmolCr	
			1-99 Years: 0.0-1.2 umol/mmolCr	
Fumarate, Urine:			0-29 Days: 1.0-26.4 umol/mmolCr	
			30-364 Days: 2.2-19.8 umol/mmolCr	
			1-99 Years: 0.3-2.6 umol/mmolCr	
Uracil, Urine:			0-29 Days: 0.0-1.6 umol/mmolCr	
			30-364 Days: 0.0-7.7 umol/mmolCr	
			1-99 Years: 0.0-5.1 umol/mmolCr	
MethylSuccinate, Ur:			0-29 Days: 0.0-5.6 umol/mmolCr	
			30-364 Days: 0.1-6.8 umol/mmolCr	
			1-99 Years: 0.0-1.4 umol/mmolCr	
Succinate, Urine:			0-29 Days: 1.1-219.5 umol/mmolCr	
			30-364 Days: 7.7-189.6 umol/mmolCr	
			1-99 Years: 0.3-27.4 umol/mmolCr	
EthylMalonate, Urine:			0-29 Days: 0.6-34.5 umol/mmolCr	
			30-364 Days: 0.7-82.7 umol/mmolCr	
			1-99 Years: 0.5-6.2 umol/mmolCr	
Benzoic Acid, Urine:			0-29 Days: 0.0-9.6 umol/mmolCr	
			30-364 Days: 0.0-16.3 umol/mmolCr	
			1-99 Years: 0.0-14.6 umol/mmolCr	
MethylMalonate, Urine:			0-29 Days: 0.0-2.0 umol/mmolCr	
			30-364 Days: 0.0-2.2 umol/mmolCr	
			1-99 Years: 0.0-0.6 umol/mmolCr	
3-HydroxyIsoValerate:			0-29 Days: 0.0-72.4 umol/mmolCr	
			30-364 Days: 1.7-119.3 umol/mmolCr	
			1-99 Years: 2.1-27.3 umol/mmolCr	
Malonate, Urine:			0-29 Days: 0.0-0.4 umol/mmolCr	
			30-364 Days: 0.0-0.4 umol/mmolCr	
			1-99 Years: 0.0-0.1 umol/mmolCr	
3OH2MethButyrate, Ur:			0-29 Days: 0.0-0.8 umol/mmolCr	
			30-364 Days: 0.0-5.6 umol/mmolCr	
			1-99 Years: 0.0-1.3 umol/mmolCr	
AcetoAcetate, Urine:			0-29 Days: 0.0-0.1 umol/mmolCr	
			30-364 Days: 0.0-2.3 umol/mmolCr	
			1-99 Years: 0.0-0.5 umol/mmolCr	
2OH-IsoValerate, Ur:			0-29 Days: 0.0-0.2 umol/mmolCr	
			30-364 Days: 0.0-0.1 umol/mmolCr	
			1-99 Years: 0.0-0.1 umol/mmolCr	
3HydroxyButyrate, Ur:			0-29 Days: 0.1-4.4 umol/mmolCr	
			30-364 Days: 0.4-9.9 umol/mmolCr	
			1-99 Years: 0.1-2.6 umol/mmolCr	
Oxalic Acid, Urine:			0-29 Days: 2.2-73.4 umol/mmolCr	
			30-364 Days: 6.9-76.4 umol/mmolCr	
			1-99 Years: 0.7-12.4 umol/mmolCr	
2HydroxyButyrate, Ur:			0-29 Days: 0.2-9.6 umol/mmolCr	
			30-364 Days: 0.0-15.0 umol/mmolCr	
			1-99 Years: 0.0-2.7 umol/mmolCr	
Lactate, Urine:			0-29 Days: 35.5-282.4 umol/mmolCr	
			30-364 Days: 15.4-198.6 umol/mmolCr	
			1-99 Years: 2.9-47.2 umol/mmolCr	

Test Changes (cont.)

Test Name	Order Code	Billing Code	Change	Effective Date
Oxycodone Confirmation, Urine	UOXYCC	82615	CPT: 80365x1	1/5/2015
Perphenazine	PRPHEN	PRPHEN	Special Information: Draw specimen prior to next dose - at steady state concentration. Please provide the following information if available 1.) Dose - List drug amount and include the units of measure 2.) Route - List the route of administration (IV, oral, etc.) 3.) Dose Frequency - Indicate how often the dose is administered (per day, per week, as needed, etc.) 4.) Type of Draw - Indicate the type of blood draw (Peak, Trough, Random, etc.) Component: Test build may need to be modified.	2/17/2015
Phencyclidine Confirmation, Urine	UPCPC	87652	CPT: 83992x1	1/1/2015
Propafenone	PROPA	76140	Special Information: Draw specimen prior to next dose - at steady state concentration. Please provide the following information if available: 1.) Dose - List drug amount and include the units of measure 2.) Route - List the route of administration (IV, oral, etc.) 3.) Dose Frequency - Indicate how often the dose is administered (per day, per week, as needed, etc.) 4.) Type of Draw - Indicate the type of blood draw (Peak, Trough, Random, etc.) Component: Test build may need to be modified.	2/17/2015
Sm/RNP Antibody	NRNP	83705	Methodology: Immunoassay (IA) Clinical Information: Smith (Sm)/U1-RNP Antibody is detected in patient with mixed connective tissue disease having features of systemic lupus erythematosus (SLE) Days Performed: Sunday-Friday Reference Range: NRNP/SM IgG Autoantibodies: < 1.0 Units Stability: Frozen: 30 Days Refrigerated: 7 Days Ambient: 4 Days Primary Name: Sm/RNP Antibody	1/6/2015
Sotalol	SOTAL	89449	Special Information: Draw specimen prior to next dose - at steady state concentration. Please provide the following information if available: 1.) Dose - List drug amount and include the units of measure 2.) Route - List the route of administration (IV, oral, etc.) 3.) Dose Frequency - Indicate how often the dose is administered (per day, per week, as needed, etc.) 4.) Type of Draw - Indicate the type of blood draw (Peak, Trough, Random, etc.)	2/17/2015
Streptozyme	STRPTO	79194	Reference Range: Streptozyme: None Detected	2/17/2015
Synthetic Cannabinoid Metabolite, Screen with Confirmation	K2	89621	CPT: 80302x1	1/1/2015
Tapentadol and Metabolite Confirm/Quantitation, Urine	TAPENU	89646	CPT: 80372x1	1/1/2015
Testosterone, Free, Adult Males by ED/LC-MS/MS	FTESAM	89350	Days Performed: Sunday, Wednesday-Saturday	2/17/2015
Testosterone, Free/Total, Males by ED/LC-MSMS	FTTESM	89286	Days Performed: Sunday, Wednesday-Saturday	2/17/2015
Tobramycin, Post Dose	TOBRPO	52022	Specimen Requirements: 1.0 mL Plasma - sodium or lithium heparin (Green - 0.5 mL minimum); Refrigerated. Collect 30 minutes after completion of infusion. Centrifuge, then transfer plasma to a clean, tightly sealed tube and refrigerate. Freeze if storage/transport time will be longer than 24 hours. Alternative - 1.0 mL Serum - SST (Gold - 0.5 mL minimum); Refrigerated. Methodology: Fluorescence Polarization Immunoassay (FPIA) Reported: 8 Hours Stability: Frozen: 1 Month Refrigerated: 3 Days Ambient: 4 Hours	4/2/2015

Test Changes (cont.)

Test Name	Order Code	Billing Code	Change	Effective Date
Tobramycin, Pre Dose	TOBRPR	52020	Specimen Requirements: 1.0 mL Plasma - sodium or lithium heparin (Green - 0.5 mL minimum); Refrigerated. Collect 5-90 minutes before next infusion. Centrifuge, then transfer plasma to a clean, tightly sealed tube and refrigerate. Freeze if storage/transport time will be longer than 24 hours. Alternate - 1.0 mL Serum - SST (Gold - 0.5 mL minimum); Refrigerated. Methodology: Fluorescence Polarization Immunoassay (FPIA) Reported: 8 Hours Stability: Frozen: 1 Month Refrigerated: 3 Days Ambient: 4 Hours	4/2/2015
Tobramycin, Random	TOBRRA	76078	Specimen Requirements: 1.0 mL Plasma - sodium or lithium heparin (Green - 0.5 mL minimum); Refrigerated. Centrifuge, then transfer plasma to a clean, tightly sealed tube and refrigerate. Freeze if storage/transport time will be longer than 24 hours. Alternate - 1.0 mL Serum - SST (Gold - 0.5 mL minimum); Refrigerated. Methodology: Fluorescence Polarization Immunoassay (FPIA) Reported: 8 Hours Stability: Frozen: 1 Month Refrigerated: 3 Days Ambient: 4 Hours	4/2/2015
VWF Exon 28 Genotyping*	VWEX28	90426	CPT: 81403x1, G0452x1	2/2/2015

*=Test Modification from January Technical Update

New Tests

Test Name	Order Code	Billing Code	Test Information	Effective Date
Anti-C1q Antibody, IgG	AC1QGG	90476	Specimen Requirements: 0.5 mL Serum - SST (Gold - 0.15 mL minimum); Refrigerated. Separate serum from cells ASAP or within 2 hours of collection. Special Information: Assess risk for lupus nephritis and global SLE disease activity. Methodology: Semi Quantitative Enzyme Linked Immunosorbent. Clinical Information: The presence of the anti-C1q IgG antibody may be associated with increased risk of lupus nephritis or with systemic lupus erythematosus (SLE) global activity. Anti-C1q antibodies are not specific for SLE; strong clinical correlation with disease is recommended. Days Performed: Monday Reported: 2-16 Days CPT: 83516x1 Price: \$110.00	3/12/2015
Reticulocyte, Hemoglobin	RTIRHB	90478	Includes: - Reticulocyte - Absolute reticulocyte - Immature reticulocyte - Reticulocyte hemoglobin equivalent Specimen Requirements: 2.5 mL Whole blood - EDTA (Lavender); Refrigerated. Methodology: Automated Cell Counter Days Performed: Sunday-Saturday Reported: 8 Hours CPT: 85046x1 Price: \$31.00	2/16/2015

Discontinued Tests

Test Name	Order Code	Billing Code	Test Information	Effective Date
Adenovirus Antibody, IgG	ADNOIG	89453	Test discontinued.	12/30/2014
Adenovirus Antibody, IgM	ADNOIM	89452	Test discontinued.	12/30/2014
Catecholamines Fractionated Tilt	TPLCAT	41000	Test temporarily discontinued.	2/10/2015
Hepatitis C Virus RNA Quantitative bDNA	HCBDNA	89281	Test discontinued. HCQPCR recommended.	2/17/2015
HIV 1 Quantitative bDNA	HIVBDN	89550	Test discontinued. HIVRNA recommended.	2/17/2015
Paraldehyde & Acetaldehyde	PARACE	75683	Test discontinued.	2/26/2015